



VINOBLE COSMETICS

SAFETY PRECAUTIONS

- Users of electronic devices such as pacemakers, electrocardiogram scanners, artificial heart-lung machines, infusion pumps, etc., should not use this device. Misuse of medical electronic devices can lead to accidents and physical discomfort.
- Do not use this product near a water-filled bathtub and/or shower tray. Keep the product and its power source in a dry state.
- Do not use the device under an oven or air conditioning heater to avoid accidents or product malfunctions.
- If the product, remote control, or power cable is damaged, do not use it to avoid the risk of electric shock.
- Do not place it in areas that are easily accessible to infants and children to prevent dangerous accidents.
- This product is not waterproof, do not immerse it in water to avoid malfunctions or accidents. If it falls into water, do not continue using it.
- This product is intended for use on the face only; do not use it for purposes other than its intended use.
- Do not disassemble the device or make any modifications or repairs to ensure the normal use of the product.
- Do not use in combination with light-sensitive ingredients, acidic skincare products, or chemical scrubs.

Please refrain from treatment in the following cases:

- Severely inflamed areas with acne, etc.
- Severely inflamed areas caused by cosmetics
- Large areas of diseased or injured skin
- Skin spots due to itching, redness, or physical irritation
- Sunburned skin
- Oculodermal melanocytosis (bluish-black discoloration of the skin in the forehead, periorbital region, and cheek area)



VINOBLE COSMETICS

The following individuals are not suitable for using this product, and excessive use may lead to accidents or skin issues:

- Individuals under 18 years old
- Individuals with heart conditions
- Individuals with malignant tumors
- Individuals with keloid scars
- Individuals with severe melasma on the face
- Individuals with mental disorders
- Individuals with acute infectious diseases
- Individuals with sunburned skin
- Individuals with abnormal blood pressure
- Pregnant and breastfeeding women
- Cancer patients, especially those with skin cancer and precancerous moles
- Individuals with impaired wound healing (e.g. diabetes)
- Individuals with autoimmune diseases
- Individuals who have had Botox injections in the last 3 months
- Individuals with contact dermatitis due to synthetic fibers
- Individuals who have been using steroid medications long-term or have capillary disorders due to liver dysfunction
- If the treatment area includes permanent or temporary subcutaneous fillers, consult a doctor before use
- If undergoing hospital treatment, consult a doctor before use

Note:

- Ensure your hands are dry when operating this product.
- This product is a precision instrument. Do not fold it, pull on it, sit on it, or cause any short circuit that could impair its use. Please store it properly after use.
- For cleaning, use disinfecting wipes.